

## Appendix G: Instructions for Use

The following steps describe the process of using the Weight Wizard:

Step-1: Turn on the scale using the button on it.

Step-2: Place the display in your desired location where it can be clearly seen. The display has a movable stand that can be mounted anywhere on a flat surface (Figure #).




Figure #: Wireless display with mountable case

Step-3: Step on the scale (Figure ##).




Figure ##: Person stepping on a bathroom scale

Step-4: Perform the activities while you are on the scale following the occupational therapist or clinical therapist instructions.

Step-5: Read the value of weight bearing on the display scale.